




# Zoll OneStep Pad Placement

# Zoll OneStep Adult Pads




# Standard Pad Placement

Appropriate pad placement for the majority of adult code events




# Anterior/Posterior Pad Placement

Can be considered for cardioversion



# Standard Placement with CPR Puck Detached

Can be used for patients with increased body habitus to ensure proper pad and puck placement.



# CPR Puck

- Puck provides real time feedback on compression quality
- Zoll will give suggestions based off the CPR puck feedback (i.e. "Push harder" or "Good Compressions")
- This feedback ensures high quality CPR is being performed in a code event as well as provides data for debriefing and quality assurance.